

Meso Skin Perfector®

EN-2-2020



**MESO SKIN PERFECTOR®**  
**Step by Step**

**mesoskinline®**

MILANO - HAMBURG - COPENHAGEN *By Aleksandra Kjaersfeldt*



# CONTOURING YOUR FACE

## MAKING MAGIC WITH MESO SKIN PERFECTOR®

### HOW TO USE

#### MESO SKIN PERFECTOR® LIGHT

Recommended for people with very light skin tones seeking the reduction of freckles and skin discolouration caused by acne and ageing.

Use this shade for a facial meso-contour BB treatment as a highlighter (see picture).



#### MESO SKIN PERFECTOR® NATURAL

Medium beige with rose undertones which is ideal for people with warmer skin tones. Use it to achieve a balanced skintone for pale skin with brown hyperpigmentation.



#### MESO SKIN PERFECTOR® MEDIUM

Ideal for tanned skin or people with bronze skin tones, and we recommend it for those with medium to dark complexions. Use this shade as a dark contour for a facial meso-contour treatment (see picture).



Recommended Meso Handpiece



Recommended Meso Microlifting Module



#### For optimal and risk-free results

- 4-8 treatments with approx. one week apart, depending on age and skin type.
- Maintained with a single treatment every 3 months.

# PROFESSIONAL USE: INSTRUCTIONS STEP BY STEP

## MAKING MAGIC WITH MESO SKIN PERFECTOR®

4-8 treatments with approx. one week apart, depending on age and skin type. Maintained with a single treatment every 3 months.

### STEP 1: MESO PURE PEEL

Prepare the skin properly with MESO PURE PEEL before undergoing MESO SKIN PERFECTOR® treatment.

- Cleanse the area to be treated, carefully removing the makeup
- Pre-procedure cleansing solution / wipe (Meso Pure Peel)

#### Active ingredients

Lactic Acid, Citric Acid, Malic Acid, Bisabolol, Citrus Grandis Fruit Extract, Vaccinum Myrtillus Leaf Extract, Citrus Lemon Peel Extract



### STEP 2: SELECT THE CORRECT NEEDLE MODULE

Preparing the handpiece with Microlifting Module or similar.

**NEVER** with Microneedling Modules.

Adjust Microlifting Module to full length and set the machine to a speed 100-120 depending on skin sensitivity.

You can treat all areas of the hands, face, neck, décolleté or body with the microlifting module, as it is non-invasive, will not traumatize the skin but will improve skin quality even on a scar tissue.



### STEP 3: APPLYING MESO SKIN PERFECTOR®

#### Select the correct color

- **Light:** for light skin tones
- **Natural:** for rose skin tones
- **Medium:** for darker skin tones

Apply a thin layer of MESO SKIN PERFECTOR® to a SMALL area at a time. Treat each area (5 cm x 10 cm) in circular movements with a moderate pressure on cheeks and no pressure around the eyes or forehead for 2-3 minutes. Treat pigment spots or scar tissue a little longer. Do not rinse the skin; leaving the residual of the MESO SKIN PERFECTOR® on the skin proceed to the next step of the treatment.



### STEP 4: MESO ANTI-AGE MASK

Post-mesotherapy transdermal patch mask to be applied directly on the skin which still has the residual of the MESO SKIN PERFECTOR® from the previous treatment step. This helps the skin absorb even more of the product. Leave the mask on for 20-40 minutes depending on clients age and skin condition, as mask has strong anti-age action and skin rejuvenation properties. Hydrogel Face Mask, easy and fast to use (pre-cut). It's 3D reticular structures ensures maximum adhesion to the skin. It is dampened with active ingredients specially selected for treatment of different types of skin. Soothing, moisturizing the skin and does not cause irritation.



### STEP 5: MESO LIFT AND PROTECT (SPF50+)

After removing the mask, wipe the rest of the MESO SKIN PERFECTOR® off with water.

Apply a thin layer of MESO LIFT AND PROTECT on the treated area. Follow up with a light mineral powder if necessary. MESO LIFT AND PROTECT is a photoprotector, suitable for everyday use. Ideal day cream for dry hyper-reactive skin, photodermatosis, photosensitivity.

MESO LIFT AND PROTECT is a multifunctional anti-age product with a strong action lasting up to 8 hours.



**Headquarter, sales, distribution and showroom**

**mesoskinline**<sup>®</sup>  
mesoskinline ApS

Juelsmindevej 57  
DK-7120 Vejle  
Denmark

(+45) 71 74 11 08  
sales@mesoskinline.com  
www.mesoskinline.dk

Innovative  
non-invasive  
mesotherapy solutions



**mesoskinline**<sup>®</sup>

MILANO - HAMBURG - COPENHAGEN *By Aleksandra Kjaersfeldt*